Pega Medical

Fassier — Duval Telescopic Im System's Rescue System 2.0



The Rescue System was conceived to retrieve components of the Fassier-Duval telescopic IM System after the treatment has been completed or in case of revision surgery.

SURGICAL TECHNIQUE

FASSIER – DUVAL TELESCOPIC IM SYSTEM'S RESCUE SYSTEM 2.0

Surgical Technique



IN ORDER TO RETRIEVE OR REVISE AN IMPLANT, NO PERMANENT DEFORMATION OF THE NAIL SHOULD HAVE OCCURRED; OTHERWISE THIS METHOD WILL NOT BE APPROPRIATE (SEE STEP 1 OF MALE RETRIEVER TECHNIQUE).

TABLE OF CONTENTS

| Sterilization and maintenance of the Rescue System | 2 |
|--|---|
| Assembly of the Male Retriever | 3 |
| Female Retrieval technique | 5 |
| Male Retrieval technique | 6 |



STERILIZATION AND MAINTENANCE OF THE RESCUE SYSTEM

The Male Retriever should be sterilized with all components disassembled as shown in the picture bellow. Sterilization should be performed by following the information included in the Guidance for Instrument Care and instructions for Use pamphlet.

This instrument should be cleaned MANUALLY before it is cleaned with ultrasound. Completely dry the instrument once cleaned by with the use of compressed air or drying oven.



Tubes

ASSEMBLY OF THE MALE RETRIEVER

Assembly of the Male Retriever components should be done after sterilization. The size of the shaft and tube are identified by the size of the implant to be retrieved. The sizes of each component of the Male Retriever are engraved on each component and are indicated by a color code as shown on the table below. Select the required size for use before starting the assembly process.

| CATALOG # | IMPLANT SIZE | COLOR CODE |
|-----------|--------------|------------|
| MR232 | 3.2 mm | Yellow |
| MR240 | 4.0 mm | Red |
| MR248 | 4.8 mm | Blue |
| MR256 | 5.6 mm | Black |
| MR264 | 6.4 mm | Rust |

STEP 1





ASSEMBLY OF THE MALE RETRIEVER

STEP 2

Verify that the flats on the driveshaft are aligned with the flats on the handle. If not, use the torque handle to turn the Driveshaft in a **counter-clockwise** direction until the flats line-up.

Insert the shaft and tube assembly into the handle. Flats should line-up.



Step 3

Slide the cap along the tube lining-up the flat surfaces. The cap will « click » into position.



Female Driver

Once the Female Retriever is engaged

turn locking nut to lock pliers

STEP 1

Use the Female Driver to unscrew the female component.

The Female Driver must be fully entered into the hexagonal drive of the implant. The instrument may break if rotated when partially entered or miss-aligned.

If the female component cannot be removed with the Female Driver use the Female Retriever to pull it out (see step 2).

If bone has grown on top of the female head, reaming with a drill or trephine will be necessary to expose the female component.

| IMPLANT SIZE | DRILL/TREPHINE |
|--------------|----------------|
| 3.2 mm | Ø 7.0 mm |
| 4.0 mm | Ø 8.5 mm |
| 4.8 mm | Ø 10.0 mm |
| 5.6 mm | Ø 10.75 mm |
| 6.4 mm | Ø 11.5 mm |

STEP 2

Insert the Female Retriever in the hexagon of the female component as shown in the picture.

Once the instrument is engaged, pull while applying a counterclockwise rotation.

STEP 1

If the implant has been bent due to a recurrence of deformity or refracture, the nail will have to be removed through an osteotomy at the apex of the bend after cutting the male component in 2 segments: proximal and distal to the bend.

Both the Male Driver and Male Retriever can be used for removal of the male component. If the Male Driver is used, it must be inserted through the intramedullary canal until it is fully engaged over the wings of the male component. It is important to note that due to patient growth, the Male Driver may be too short to reach the implant wings.

Partial engagement of the Male Driver onto the male component wings may lead to breakage of the instrument due to increased stress on the instrument. If proper engagement of the Male Driver cannot be achieved, it is recommended to use the Male Retriever to remove the male component (see step 2).

Male Driver

Step 2

The Male Retriever has a larger outer diameter than the female component. Therefore, once the female has been retrieved, reaming with the appropriate drill size to the tip of the male nail may be necessary before the tool can be inserted into the same canal. The table below shows the reamer sizes required for each size of implant to be removed.

| CATALOG # | IMPLANT SIZE |
|-----------|--------------|
| MR232 | DR148 |
| MR240 | DR148 |
| MR248 | DR156 |
| MR256 | DR164 |
| MR264 | DR164 |



Β

STEP 3

Under fluoroscopy, the end of the male component should be located and the tip of the Male Retriever tool inserted over it. **Make sure that the end of the male component is securely inserted in the hole of the Male Retriever shaft before tightening.**

The hole in the tube should not be visible.





Tightening of the Male Retriever handle without the male component inside will result in damage to the shaft and render it unusable. A new Shaft will be necessary since they are single use instruments.



To lock the shaft over the male component, turn the torque handle **counterclockwise** until the torque limit is reached. Remove the torque handle after tightening.

B Once the nail is secured in the tool, turn the Male Retriever handle in a counterclockwise direction. Simultaneous rotation and pulling might be necessary as the male component is unscrewed.

C To remove the retrieved implant from the Rescue Tool, press the Release Button all the way down using the Release Key.





MALE RETRIEVAL TECHNIQUE

STEP 4

Disassembly of the shaft and tube from the handle is done by pulling the components upwards after the cap has been removed.



While maintaining the Release Key engaged into the Release Button rotate the driveshaft clockwise until the flats line-up. The shaft and tube assembly can then be pulled upwards.





Pega Medical

1111 Autoroute Chomedey, Laval, Quebec CANADA H7W 5J8 Phone: 450-688-5144 • Fax: 450-688-1977 info@pegamedical.com Pegamedical.com

© 2015 Pega Medical, Inc.

Distributed by